



## **CASE STUDY:** **FROM NEAR-AMPUTATION TO RECOVERY** **Puracyn<sup>®</sup> Plus Professional Formula in Tonga**

**WOUND IRRIGATION**  
Non-sensitizing, Non-toxic, Non-irritating

### **Dedicated to a Radical Change**

In partnership with 'Amanaki Fo'ou, a Utah-based team of diabetes specialists, Michelle Moore, MSN, RN, CWCN, WCC, CLWCP, and her Puracyn<sup>®</sup> Team members traveled to Tonga to support wound care for people suffering from diabetes and chronic wounds.

### **Case Study**

A 37-year-old man with a history of venous insufficiency and obesity presented with a right leg covered with necrotic tissue, exposed tendon, and heavy exudate. Facility clinicians requested amputation.



#### **Patient at Initial Presentation, 6/26**

After carefully examining the leg, the Puracyn Team suggested the following protocol in lieu of amputation: Surgical debridement, followed by a daily regimen of soaking the wound in Puracyn Plus for ten minutes, dressing the wound in Puracyn Plus-soaked gauze, and wrapping with compression bandages.



#### **Wound Progression**

The Puracyn Team offered nutrition and health counseling during the course of treatment. As a result, over the two-and-a-half-month treatment period, the patient successfully lost 100 pounds.



#### **Wound Progressing Toward Closure, 8/29**

After seven weeks of daily soaking and wrapping, the patient was discharged. The protocol changed to soaking with Puracyn Plus and application of an Unnas Boot with dressing changes once every week at the clinic until closure.



#### **Conclusion**

Daily soaking with Puracyn Plus, in conjunction with appropriate debridement and dressings, help non-healing wounds achieve closure.

**VISIT [WWW.PURACYNPRO.COM](http://WWW.PURACYNPRO.COM) TO LEARN MORE**



Clear the way.

## Expanding the Possibilities of Wound Care

Hypochlorous technology is both radically simple and based on powerful science. It's the foundation of every Puracyn® Plus Professional Formula product and the best strategy when it comes to wound care protocols. Safe, effective, and easy to use, our suite of products promote healing, leading to improved clinical and economic outcomes.

Learn more about the Puracyn® Plus Professional Formula Family of Products at [puracynpro.com](http://puracynpro.com) or call 866.318.3116 to talk to a representative today.



**Puracyn**  
plus

PROFESSIONAL FORMULA